

	<b>Cincinnati Fire Department Fire Training Supplement DRILL BOOK</b>	<b>SECTION #6 EMS Equipment and Procedures</b>
<b>Date:</b> January 2006	<b>TOPIC TITLE:</b> SurgiLance Lancet	<b>Total Pages: 1</b>
<b>Section #: 6</b>		<b>Topic #: 2</b>

## TOPIC #2: SurgiLance®



### Description:

One-step single-use safety lancet for blood sampling

### Features:

1. The “one-step” needle is concealed before and after use. It is automatically rendered inoperative after use.
2. The “one-step” safety lancet is pre-armed and ready for use.
3. The lancet is triggered by a gentle push on the test site area.

### Use:

1. Clean test site with alcohol pads and wipe “DRY” with a sterile 4x4 gauze pad to prevent the alcohol on the skin from affecting the reading obtained by the Glucometer.
2. Hold the safety lancet in your hand. (picture #1)



3. Apply gentle pressure on the test site area, you will hear a click sound indicating the needle has fired. (picture #2)
4. Remove the safety lancet from the skin and dispose of immediately into an approved sharps disposal container.
  - a. Sharps Shuttle – SINGLE USE ONLY
  - b. Flip-top sharps container
  - c. 1 quart sharps container
5. Obtain blood sample with the Glucometer.

